





Quickarrays, Inc

Phone: 415-863-2380
Support: info@quickarrays.com
Web: www.quickarrays.com

Advancing Biomedical Science Through Tissue Arrays

Cat No: **BRC711** – Breast cancer tissue array

Lot#	Cores	Size	Cut	Format	QA/QC	
BRC71101	71	1.1mm	4um	8X9	H&E, IHC anti- Cytokeratin	 

Recommended applications: For Research use only. RNA or protein breast cancer/non-tumor tissue profiling using IHC or ISH; Antibody characterization.

Description: Breast cancer tissue array, 71 cores including 71cases of breast cancers (some cases may run out tumor cells due to small core size).

All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They can be stored for use at 4C for up to six months from the date of shipment. **There may be 5 to 10% of tissue core loss.**

Array Position	Sex	Age	Anatomic Site	Pathology	Grade	Stage
A01	F	37	Breast	Invasive ductal carcinoma, few tumor cells	I	T3N1M0
A02	F	54	Breast	Invasive ductal carcinoma	II	T2N0M0
A03	F	43	Breast	Invasive ductal carcinoma	II	T3N1M0
A04	F	38	Breast	Invasive ductal carcinoma	II	T1N0M0
A05	F	58	Breast	Fat tissue	null	null
A06	F	54	Breast	Invasive ductal carcinoma	III	T2N1M0
A07	F	58	Breast	Invasive ductal carcinoma	II~III	T3N1M0
A08	F	66	Breast	Invasive ductal carcinoma	II	T3N1M0
A09	F	40	Breast	Invasive ductal carcinoma	III	T1N0M0
B01	F	40	Breast	Invasive ductal carcinoma	II	T2N1M0
B02	F	29	Breast	Invasive ductal carcinoma	II~III	T2N1M0
B03	F	46	Breast	Invasive ductal carcinoma	II	T2N0M0
B04	F	39	Breast	Invasive ductal carcinoma	III	T2N1M0
B05	F	40	Breast	Invasive ductal carcinoma	II	T3N0M0
B06	F	44	Breast	Invasive ductal carcinoma	II~III	T2N0M0
B07	F	46	Breast	Invasive ductal carcinoma	II~III	T2N0M0
B08	F	35	Breast	Invasive ductal carcinoma	I~II	T2N1M0
B09	F	55	Breast	Invasive ductal carcinoma	II	T2N1M0
C01	F	38	Breast	Invasive ductal carcinoma	II	T3N1M0
C02	F	35	Breast	Invasive ductal carcinoma	II~III	T2N0M0
C03	F	40	Breast	Invasive ductal carcinoma	II~III	T2N0M0
C04	F	50	Breast	Mucinous adenocarcinoma	II~III	T2N0M0
C05	F	40	Breast	Invasive ductal carcinoma	II	T2N0M0
C06	F	53	Breast	Ductal carcinoma in situ	I	TisN0M0
C07	F	49	Breast	Invasive ductal carcinoma	II	T3N1M0
C08	F	45	Breast	Ductal carcinoma in situ	I	TisN0M0
C09	F	62	Breast	Invasive ductal carcinoma	III	T1N0M0
D01	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
D02	F	33	Breast	Invasive ductal carcinoma	II	T2N0M0
D03	F	37	Breast	Invasive ductal carcinoma, few tumor cells	II~III	T3N1M0
D04	F	47	Breast	Invasive ductal carcinoma	III	T1N0M0
D05	F	50	Breast	Invasive ductal carcinoma	II~III	T1N0M0
D06	F	64	Breast	Invasive tubular carcinoma	II~III	T3N1M0
D07	F	44	Breast	Invasive ductal carcinoma	II	T3N1M0
D08	F	55	Breast	Invasive ductal carcinoma	II	T2N1M0
D09	F	47	Breast	Invasive ductal carcinoma	II~III	T2N0M0
E01	F	50	Breast	Invasive ductal carcinoma	II~III	T2N1M0



Quickarrays, Inc

Phone: 415-863-2380
Support: info@quickarrays.com
Web: www.quickarrays.com

Advancing Biomedical Science Through Tissue Arrays

E02	F	51	Breast	Invasive ductal carcinoma	II	T4N1M0
E03	F	51	Breast	Invasive ductal carcinoma	II~III	T3N1M0
E04	F	36	Breast	Invasive ductal carcinoma	II~III	T3N1M0
E05	F	36	Breast	Invasive ductal carcinoma	II~III	T2N0M0
E06	F	43	Breast	Invasive ductal carcinoma	II	T4N1M0
E07	F	50	Breast	Invasive papillary carcinoma	II~III	T1N0M0
E08	F	47	Breast	Invasive ductal carcinoma	I~II	T3N1M0
E09	F	48	Breast	Invasive ductal carcinoma	II	T2N0M0
F01	F	66	Breast	Invasive ductal carcinoma, tubular type	I~II	T2N0M0
F02	F	50	Breast	Invasive ductal carcinoma	II	T2N0M0
F03	F	66	Breast	Invasive ductal carcinoma	III	T2N0M0
F04	F	49	Breast	Invasive ductal carcinoma	II	T2N0M0
F05	M	61	Breast	Ductal carcinoma in situ, papillary	I~II	TisN0M0
F06	F	33	Breast	Invasive ductal carcinoma	III	T2N0M0
F07	F	78	Breast	Invasive ductal carcinoma	II~III	T4N0M0
F08	F	40	Breast	Invasive lobular carcinoma	II~III	T2N0M0
F09	F	68	Breast	Invasive ductal carcinoma	II~III	T2N0M0
G01	M	82	Breast	Invasive ductal carcinoma, few tumor cells	II~III	T2N0M0
G02	F	75	Breast	Invasive ductal carcinoma	II~III	T2N0M0
G03	F	48	Breast	Invasive ductal carcinoma	II~III	T4N1M0
G04	F	77	Breast	Invasive ductal carcinoma	III	T2N0M0
G05	F	46	Breast	Invasive ductal carcinoma	II~III	T2N0M0
G06	F	53	Breast	Normal breast	null	null
G07	F	43	Breast	Ductal carcinoma in situ	I	TisN0M0
G08	F	61	Breast	Normal breast	null	null
H01	F	40	Breast	Invasive ductal carcinoma	I	T2N0M0
H02	F	79	Breast	Invasive ductal carcinoma	II	T2N1M0
H03	F	55	Breast	Invasive ductal carcinoma	II	T3N1M0
H04	F	51	Breast	Invasive lobular carcinoma	II~III	T2N0M0
H05	F	54	Breast	Invasive ductal carcinoma	III	T2N0M0
H06	F	43	Breast	Invasive lobular carcinoma	I	T4N1M0
H07	F	57	Breast	Invasive ductal carcinoma	II~III	T2N0M0
H08	F	47	Breast	Invasive ductal carcinoma	I~II	T2N0M0
H09	F	32	Breast	Invasive ductal carcinoma	I~II	T2N0M0

Notes: Bake at 60C for ~60 minutes before use. If antigen retrieving is needed, it is important to avoid **direct-boiling and high pH or high strength** antigen retrieving buffer.

Certified by: Langxing Pan, M.D.

TNM Classification: Breast carcinoma

T- Primary tumor

TX - Primary tumor cannot be assessed

TO - No evidence of primary tumor

Tis - Carcinoma in situ, intraductal carcinoma or lobular carcinoma in situ, or Paget's disease of the nipple with no tumor

T1 - Tumor 2 cm or less in greatest dimension

T2 - Tumor more than 2 cm but not more than 5 cm in greatest dimension

T3 - Tumor more than 5 cm in greatest dimension

T4 - Tumor of any size with direct extension to chest wall or skin



Quickarrays, Inc

Phone: 415-863-2380

Support: info@quickarrays.com

Web: www.quickarrays.com

N - Regional lymph nodes

NX - Regional lymph nodes cannot be assessed

NO - No regional lymph node metastasis

N1 - Metastasis to movable ipsilateral lymph node(s)

N2 - Metastasis to ipsilateral lymph nodes(s) fixed to one another or to other structure

N3 - Metastasis to ipsilateral internal mammary lymph node(s)

M - Distant metastasis

MX - Distant metastasis cannot be assessed

MO - No distant metastasis

M1 - Distant metastasis

TNM classification of malignant tumours, Fifth Edition (1997)